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
Preventing peri-implant diseases:
the new challenge in implantology

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
Claiming «zero peri-implantitis» is an illusion. It is as if a toothpaste could guarantee «zero cavities». However, a therapeutic procedure aimed at restoring smile and chewing function can sometimes, paradoxically, be the cause of a new pathology.

Today, discussions on implant osseointegration also emphasize that bone stability is closely linked to soft tissue health. Indeed, maintaining healthy gums is more important than prioritizing excessive aesthetics when designing a crown.

Thanks to pioneers like Prof. Sami SANDHAUS and Dr. André BENHAMOU, an intuition was born in the 1980s: to bring ceramics into contact with soft tissue in implantology. The European Association for Osseointegration (EAO), in publishing the 2030 trends in modern implantology (Delphi study), supports the use of biomaterials adapted to each tissue.

Prevention in implantology appears to be primarily focused on the critical area of biological width. Drawing inspiration from natural teeth and their defense mechanisms, exploring biomaterials such as Zirconia, and adapting design and colorimetry are all tools at our disposal to support functional connective tissue and the epithelial attachment layer, which are essential for peri-implant health.

This white paper is an invitation to go beyond osseointegration. Mucositis and peri-implantitis can no longer be overlooked if we are to maintain the high long-term success rates of the implant restorations we provide for our patients.

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
Prevalence of peri-implant diseases

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Peri-implant complications follow a well-established clinical pattern, with alarming incidence rates.

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Mucositis is the implant equivalent of gingivitis around a tooth. The aetiology of this inflammatory condition is inadequate oral hygiene, combined with the nature of the biomaterials in contact with the gum. Just as periodontitis follows uncontrolled gingivitis, peri-implantitis occurs when mucositis is not intercepted.

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Mucositis, which affects almost 50% of implants, is the first warning sign before an inflammatory reaction occurs in all peri-implant tissues.

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Peri-implantitis is the implant equivalent of periodontal disease on natural teeth: an aggressive inflammatory and infectious condition that causes progressive loss of the implant's supporting tissues.

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Loss of osseointegration marks the ultimate stage of implant complications: mobility and pain make implant removal inevitable, equivalent to the avulsion of a natural tooth affected by severe periodontitis.

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Importance of the epithelial-connective attachment

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The balance of the epithelial-connective tissues around natural teeth

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
Around a natural tooth, tissue integrity and stability are guaranteed by :

- **At the epithelial level** = a keratinized oral mucosa, a keratinized oral epithelium at the top of the sulcus, a non-keratinized sulcular epithelium followed by a short junctional epithelium in contact with the enamel. The hemidesmosomes ensure the adhesion of this epithelium to the tooth enamel.
- **At the connective level** = a sleeve composed of a dense network of collagen fibers that adhere to the cementum.

The primary function of this complex architecture is to seal hermetically the supracrestal zone in order to protect and preserve the bone tissue from inflammatory and infectious threats originating from the oral cavity.

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


FOCUS BIOLOGICAL WIDTH

This critical area extends from the deepest point of the sulcus to the crestal bone. A dense network of oblique and circular fibers straps all the supporting cervical tissues.

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
At the supracrestal level, the oblique fibers, anchored perpendicularly into the cement, provide a true biological attachment by preserving the periodontal health. The mechanical action of the circular fibers seals and protects the underlying periodontal tissues.

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Behaviour of connective fibers in
contact with a natural tooth

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
The combination of oblique and circular fibers acts to seal the cervical contour, maintains a functional volume of connective tissue and, most importantly, prevents epithelial invagination. The short junctional epithelium thus formed is maintained at a healthy distance from the bone crest.

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The biological width surrounding an implant results from a physical adaptation rather than an actual biological attachment. The collagen fibers are not anchored to either the implant or the prosthetic component, and are oriented randomly. Fewer circular fibers are embedded in a less dense connective tissue. A downward invagination of the junctional epithelium is observed after implant placement. This long junctional epithelium does not provide a tight barrier around the implant or the prosthetic surface.

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Behaviour of connective fibers in
contact with an implant

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Since there is no place in the body where epithelial tissue is in direct contact with bone, the observed crestal bone loss is a homothetic replica of this epithelial invagination aimed at recreating a biological width.

The weakening of the biological, physiological and mechanical defenses of this subcritical contour could at least trigger mucositis.

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

The choice of the biomaterial influences both osseointegration of the implant and architecture of the epithelial-connective attachment, essential for the stability and protection of peri-implant tissues.

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



Various supracrestal emergence geometries are now available aiming to preserve soft-tissue health.

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
SWITCHING PLATFORM PLATFORM SWITCHING

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

The translucency of titanium in transgingival applications remains a concern, and implantology manufacturers have developed industrial processes to overcome these drawbacks.

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



BIOMIMETIC
inspired by natural teeth

ANTHOGYR 
A Straumann Group Brand

BIOHORIZONS® 

 **BIOTECH DENTAL**

camlog

 **Dentsply
Sirona**

 **MEGA¹GEN**

MIS®

 **NEODENT**®
A Straumann Group Brand

 **Nobel
Biocare™**

OSSTEM®
IMPLANT

 **straumann**

 **sweden & martina**

 **TBR**®

 **TRI**
DENTAL IMPLANTS

 **ZimVie**

5 Analyse comparative des principaux designs implantaires

Comparative analysis of the main implant designs

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The interface between the dental implant and the peri-implant soft tissues plays a central role in long-term biological stability.

While most current implant systems are based on historical designs focused on osseointegration, biomimetic criteria provide a new framework based on the observation of the natural behavior of biological tissues.

Here, implants representing more than 75% of the global market* are compared according to four key criteria influencing tissue response and aesthetic integration:

- type of material in contact with bone,
- type of material in contact with soft tissue,
- supracrestal design,
- colour of the transgingival collar.


This comparison helps distinguish conventional approaches from biomimetic designs.

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**The trademarks mentioned in this document are the property of their respective owners. They are cited for illustrative purposes as part of a comparative analysis based on public technical data and scientific publications available at the time of writing.*

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
CATEGORY 2: titanium tissue level implant

Biomatériau osseux Bone biomaterial  Titanium	Biomatériau gingival Gingival biomaterial  Titanium	Bioforme Design  Tissue Level	Colorimétrie Colorimetry  Titanium
---	---	---	--

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

Trends in dental implantology up to 2030

Delphi Study - *Horizon 2030*

Identifying and predicting future trends
in implant dentistry in Europe

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The Delphi study «Horizon 2030» carried out by the European Association for Osseointegration (EAO) aims to identify and predict future trends in implant dentistry in Europe.

The overall purpose of this study is to provide information on relevant aspects of implant therapy that remain uncertain and cannot be clarified by scientific evidence.

The Delphi method involves a structured group of experts, systematic communication, individual feedback and group judgement:

- A 60-question questionnaire was sent to 138 experts selected according to geographical and professional criteria.
- Responses were collected in two rounds with feedback, then discussed at a final meeting for those points where there was no consensus (threshold set at 65%).

La prévalence de la péri-implantite va :
The prevalence of peri-implantitis will:

Augmenter
Increase


75 %

Etre similaire
Be similar

23 %

Diminuer
Decrease

2 %

Concernant le matériau de fabrication des implants dentaires, quelle sera la tendance ?
In regard to the material dental implants are made of, what do you think the trend will be?

Titane et Céramique
Titanium and Ceramic


69 %

Titane
Titanium

29 %

Céramique
Ceramic

2 %

Concernant la microtopographie de la surface des piliers, quelle sera la tendance ?
In regard to the microtopography of the abutment surface, what do you think the trend will be?

Lisse
Polished


77 %

Lisse et rugueuse
Polished and textured

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
Implantology has entered a new era in which simply restoring lost teeth is no longer sufficient. Tissue stability and the prevention of peri-implant diseases have become the true keys to the long-term success of our restorations.

Today, the architecture and physiology of the epithelial-connective tissue around natural teeth are well understood, and implantology, through concepts of varying degrees of validation, seeks to reproduce this protective barrier function.

In fact, nearly one in two implants develops mucositis, the first step towards severe complications.

The Z1 implant incorporates the following biological principles: a tissue-level design, a supracrestal zirconia collar, an anatomical emergence profile, and a natural tooth-like colorimetry, all key factors enabling biomimetic and functional epithelial-connective tissue remodeling.

By closing the door to inflammatory and microbial threats, this hybrid zirconia-titanium implant paves the way for a truly preventive approach to implantology, in which every technological choice contributes to preserving long-term tissue health.

8 Bibliographie

Bibliography

G. Da Gama Ramos, D. L. Ciotti, S. R. W. Ferreira, M. R. W. F. Margarido, R. A. Dantas, and M. N. Guerreiro, "The use of a hybrid pillar and its importance for aesthetic rehabilitation and tissue stability: a clinical report," *Case Reports in Dentistry*, vol. 2018, pp. 1–6, Jul. 2018, doi: 10.1155/2018/6850720.

M. Katafuchi, B. F. Weinstein, B. G. Leroux, Y. Chen, and D. M. Daubert, "Restoration contour is a risk indicator for peri implantitis: A cross sectional radiographic analysis," *Journal of Clinical Periodontology*, vol. 45, no. 2, pp. 225–232, Oct. 2017, doi: 10.1111/jcpe.12829.

Y.-H. Chou et al., "Association between prosthesis contour and peri-implantitis in patients compliant with supportive periodontal therapy: A retrospective cohort study," *Journal of Prosthodontic Research*, Jan. 2024, doi: 10.2186/jpr.jpr_d_23_00236.

Y. Yi, K. Koo, F. Schwarz, H. B. Amara, and S. Heo, "Association of prosthetic features and peri implantitis: A cross sectional study," *Journal of Clinical Periodontology*, vol. 47, no. 3, pp. 392–403, Jan. 2020, doi: 10.1111/jcpe.13251.

N. Mattheos, M. Janda, A. Acharya, S. Pekarski, and C. Larsson, "Impact of design elements of the implant supracrestal complex (ISC) on the risk of peri implant mucositis and peri implantitis: A critical review," *Clinical Oral Implants Research*, vol. 32, no. S21, pp. 181–202, Oct. 2021, doi: 10.1111/clr.13823.

S. Chu et al., "Restorative Emergence Profile for Single-Tooth Implants in Healthy Periodontal Patients: Clinical Guidelines and Decision-Making Strategies," *The International Journal of Periodontics & Restorative Dentistry*, vol. 40, no. 1, pp. 19–29, Dec. 2019, doi: 10.11607/prd.3697.

Molina, I. Sanz Sánchez, C. Martín, J. Blanco, and M. Sanz, "The effect of one time abutment placement on interproximal bone levels and peri implant soft tissues: a prospective randomized clinical trial," *Clinical Oral Implants Research*, vol. 28, no. 4, pp. 443–452, Mar. 2016, doi: 10.1111/clr.12818.

J. Calatrava, I. Sanz Sánchez, A. Molina, J. Bollain, C. Martín, and M. Sanz, "Effect of one time placement of the definitive abutment versus multiple healing abutment disconnections and reconnections during the prosthetic phase on radiographic and clinical outcomes: A 12 month randomized clinical trial," *Clinical Implant Dentistry and Related Research*, vol. 26, no. 5, pp. 998–1011, Jul. 2024, doi: 10.1111/cid.13361.

Su H, Gonzalez-Martin O, Weisgold A, Lee E. Considerations of implant abutment and crown contour: critical contour and subcritical contour. *Int J Periodontics Restorative Dent*. 2010 Aug;30(4):335-43. PMID: 20664835.

N. R. dos Reis et al., 'The prevalence of peri-implant mucositis and peri-implantitis based on the world workshop criteria: A systematic review and meta-analysis', *Journal of Dentistry*, vol. 160, p. 105914, Sept. 2025, doi: 10.1016/j.jdent.2025.105914.

L. J. A. Heitz-Mayfield, 'Peri-implant mucositis and peri-implantitis: key features and differences', *Br Dent J*, vol. 236, no. 10, pp. 791–794, May 2024, doi: 10.1038/s41415-024-7402-z.

T. M. A. Al-Faraj et al., "Comparison of Peri-Implant soft tissue around zirconia and titanium abutments in the aesthetic zone: A Narrative review," *Cureus*, Jul. 2024, doi: 10.7759/cureus.65782.

S. Al-Asaadi, N. Austin, P. J. Watson, D. J. Wood, A. Altaie, and F. P. Rodrigues, "Titanium-Zirconia abutment-implant assemblies: Are they alternatives for single material implants?," *Dental Materials*, Jun. 2025, doi: 10.1016/j.dental.2025.05.007.

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T. Steinherr, 'Implant Materials in Focus: A Comprehensive Review of Zirconia vs. Titanium and Their Impact on Peri-Implant Health', *JOMDR*, vol. 5, no. 4, Oct. 2024, doi: 10.52793/jomdr.2024.5(4)-79.

Kumar, S., Singha, A., Dhingra, J., Verma, G., Singh, S. D., & Kumar, A. (2024). Evaluation of zirconia and titanium abutments on peri-implant tissues: A comparative study. *Bioinformation*, 20(12), 1809–1812. <https://doi.org/10.6026/9732063002001809>

A. E. Bianchi, M. Bosetti, G. D. Jr, M. T. Sberna, F. Sanfilippo, and M. Cannas, "In vitro and in vivo follow-up of titanium transmucosal implants with a zirconia collar," *J. Appl. Biomater. Biomech.*, vol. 2, no. 3, pp. 143–50, 2004.

D.-J. Lee, J.-S. Ryu, M. Shimono, K.-W. Lee, J.-M. Lee, and H.-S. Jung, "Differential Healing Patterns of Mucosal Seal on Zirconia and Titanium Implant," *Front. Physiol.*, vol. 10, no. 796, 2019, doi: 10.3389/fphys.2019.00796.

M. Sanz et al., "European Association for Osseointegration Delphi study on the trends in Implant Dentistry in Europe for the year 2030," *Clinical Oral Implants Research*, vol. 30, no. 5, pp. 476–486, Apr. 2019, doi: 10.1111/clr.13431.

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How to recover naturally sealed tissues
in implantology?

Document non contractuel - Version 10-2025
Non contractual document

📍 24 impasse René Couzinet, 31500 Toulouse - FRANCE
☎ +33 (0)5 62 16 71 00 ✉ contact@tbr.dental
🌐 the-z1.com